

## Model 177 Irrigation Syringe

HCPCS:



- The Model 177 Irrigation Syringe is efficiently used in various E.N.T irrigation procedures
- Small, malleable tube permits deep, yet painless introduction into external ear canal
- Regulate flow by placing finger over small hole in head of instrument if air pressure is 20 psi or under, then stream will stop when finger is removed
- Used by E.N.T doctors, medical doctors, research laboratories, veterinarians
- Can be operated by the No 630 or a No 1 bulb

PRODUCTS	Description
177	Irrigation Syringe

  

ACCESSORIES	Description
177-NKHR-605	STRT FLW TIP ASSY RHODIUM PLATED 5 PK, EA
115-QB	1-OZ CLEAR GLASS, EA
177-OT-605	OUTLET TUBE ASSEMBLY 5-PK, PKG
177-QB	6-OZ CLEAR GLASS BOTTLE, EA
114-S-605	WASHER HERCULES 10-481 5-PK, PKG
630-605	CUTOFF WITH TUBING 5-PK, PKG