# **Ostomy Accessories**



# **Ostomy Accessories**

McKesson's ostomy accessories are designed for patients who have ileostomy, colostomy and urostomy stomas and need skin protection, skin cleaning, leak prevention and odor control. Our product line allows patients to continue living healthy and productive lives.

### **Features and Benefits**

### **Skin Barrier Strips**

- Keeps the pouch secure and protects the peristomal skin
- Water-resistant surface with gentle, absorbent hydrocolloid adhesive soaks up exudate and allows for use when swimming, bathing, etc.
- Pre-cut, pre-shaped strip with high-tack adhesion for easy application
- "Curved strip" available in 1/3-arc and 1/2-arc shapes to accommodate varying applications and pouching systems
- Thin-edge, no roll-up design enhances comfort and helps the strips stay smooth against the skin to allow for freedom of movement
- Ideal for patients with tape allergies or irritation
- Easy and gentle removal protects the skin from trauma caused by repetitive removal of tapes/adhesives
- Not made with natural rubber latex



### Satisfaction Guaranteed

If you are not completely satisfied with any McKesson Brands product, you may return it for a full refund or credit.

### Never Compromise on Care<sup>SM</sup>

With McKesson Brands, you'll never have to compromise. We deliver quality products at competitive prices, backed by our 100% satisfaction guarantee. We have decades of experience in the market, and we offer you support and work to help ensure you have the products you need, when you need them—making it even easier for you to provide the best care for your patients.

The product information contained in this document, including the product images and additional product materials, was collected from various supplier sources. All product claims and specifications are those of the product suppliers and have not been independently verified by McKesson Medical-Surgical or its affiliates ("McKesson"). McKesson is not responsible for errors or omissions in the product information.

The properties of a product may change or be inaccurate following the posting or printing of the product information in the document, either in the print or online version. Caution should be exercised when using or purchasing any products from McKesson's online or print documents by closely examining the product packaging and the labeling prior to use.

Due to product changes, information listed in this document is subject to change without notice. This information is placed solely for your convenience in ordering and McKesson disclaims all responsibility for its completeness and accuracy, whether or not the inaccuracy or incompleteness is due to fault or error by McKesson.

© 2018 McKesson Medical-Surgical Inc.

### **Features and Benefits**

### **Skin Barrier Rings**

- Creates a secure seal between the stoma and the ostomy appliance
- Durable design provides a secure fit, extended wear time and lasting protection against leakage
- Absorbent hydrocolloid ring remains stable when exposed to effluent
- · Excellent adhesive tack
- Conforms to create a secure "turtleneck" seal at the base of stoma
- · Moldable form requires zero cutting
- · Ideal for all types of ostomy
- Not made with natural rubber latex



### **Skin Barrier Paste**

- · Pectin-based formula
- Hydrocolloid formulation creates an absorbent, protective skin barrier
- Used as a caulk to help reduce leakage by filling creases and folds to fill in uneven skin surfaces, yielding a flatter surface for a secure fit
- Helps prevent drainage from breaching the ostomy skin barrier
- Promotes extended pouching system wear time
- Quick and reliable; easy-to-apply paste hardens upon application
- Easy-to-squeeze tube with screw-top lid
- · Not made with natural rubber latex



### **Features and Benefits**

### **Adhesive Remover**

- Helps to remove adhesive tape, adhesive residue and pouches from skin
- Cleanses skin and removes adhesives and debris from peristomal skin
- Gentle alcohol-based formula
- Packaged in convenient individually wrapped packets
- Not made with natural rubber latex



### No-Sting Skin Barrier Film

- Forms a barrier to protect the skin from adhesives, friction and body fluids
- Alcohol-free formula for sting-free comfort
- Transparent with no offensive odor
- · Increases adhesion of tapes and wafers
- Gentle formula helps prevent breakdown of normal/intact skin
- Convenient, sterile individually wrapped packets
- · Used to initiate "crusting technique"
- Not made with natural rubber latex



### **Features and Benefits**

### Skin Barrier Powder

- Absorbs moisture and enhances comfort from broken skin around the stoma or under the ostomy pouch
- Pectin-based powder forms a barrier for protection from itching, moist, irritated or broken skin
- Use with skin barrier film to complete the "crusting technique" prior to pouch application
- Helps secure tape or adhesives
- Not made with natural rubber latex

# MSKESSON Skin Barrier Powder Forms a barrier to help protect and maintain dry, adherent skin. NET WEIGHT 1 0Z DB gj MFR # 137-5724

### **Ostomy Pouch Deodorizer**

- Helps eliminate odor in ostomy pouches
- Delivers long-lasting odor protection
- · Easy to use
- No-scent formula neutralizes odor
- Available in dropper-style 1 fl. oz. bottle for added convenience
- Ideal for use with all pouching systems
- Not made with natural rubber latex



### **Features and Benefits**

### **Ostomy Lubricating Deodorant**

- Helps eliminate odor in ostomy pouches
- Lubricating surfactant makes pouch interior slippery to prevent "pancaking" and facilitates the emptying of contents
- · Clear formula reduces mess and staining
- Easy to use and provides lasting results
- Ideal for use with all pouching systems
- Available in individual packets for convenient use
- Not made with natural rubber latex



### ORDERING INFORMATION

MFR #	DESCRIPTION	QTY
176-5725	Skin Barrier Strips, 1/2 ARC, Water Resistant, Hydrocolloid, 1 in width	1 ea, 20 strips/tray, 14 trays/cs
176-5726	Skin Barrier Strips, 1/3 ARC, Water Resistant, Hydrocolloid, 1 in width	1 ea, 30 strips/tray, 14 trays/cs
176-5727	Skin Barrier Rings, 1/16 in Thickness, Absorbent, Conforming Seal, Moldable, Hydrocolloid, 2 in	1 ea, 20 rings/bx, 8 bx/cs
176-5728	Skin Barrier Film, No-sting Wipe, Sterile, 2.4 in x 2.4 in	1 ea, 25 packs/bx, 100 bx/cs
176-5729	Adhesive Remover, Wipe, Alcohol, Sterile, 2.4 in x 2.4 in	50 packs/bx, 50 bx/cs
137-5719	Ostomy Pouch Deodorizer, 1 fl oz	1 ea, 12 ea/cs
137-5720	Ostomy Pouch Deodorizer, 8 fl oz	1 ea, 6 ea/cs
137-5721	Lubricating Deodorant, 8 fl oz	1 ea, 6 ea/cs
137-5722	Lubricating Deodorant, 8 mL	1 ea, 8 ea/cs
137-5723	Skin Barrier Paste, Pectin-Based, Water Resistant, Hydrocolloid, 2 oz	1 ea, 1 ea/bx, 24 bx/cs
137-5724	Skin Barrier Powder, 1 oz	1 ea, 50 ea/cs

McKesson Medical-Surgical 9954 Mayland Drive Richmond, VA 23233 mms.mckesson.com/mckesson-brands

Post-Acute Care and Home Health 800.347.2456

Home Medical Equipment 888.822.8111

Primary Care 866.625.2679

McKesson Brands Support 800.777.4908