

Bubble Humidifier

McKesson Bubble Humidifier is used to supply humidified flow to patients undergoing oxygen therapy. Humidified oxygen flow is supplied to the patient to minimize discomfort and the potential health risks of breathing a dry flow of oxygen. Our humidifiers interface easily with all oxygen delivery equipment and wall outlet applications. The transparent bottle with water line indicators allows for easy monitoring of water levels.

Features and Benefits

- High capacity humidifier jar allows longer continuous oxygen therapy
- Maximum and minimum water level markings are clearly defined
- Pressure relief valve. Relieves pressure and sounds audible alarm when line pressure exceeds 6 psi
- Interfaces with Oxygen Concentrators, Liquid O2 Systems, O2 Tanks and Wall Outlet Applications
- Homecare, Hospital and Alternate Care Applications
- Recommended Flow Rate up to 6 LPM
- Not made with natural rubber latex



*The alarm will sound, either continuously or intermittently, when oxygen flow downstream is obstructed, and the gas pressure reaches ranges greater than 6 pounds per square inch (psi). Alarms at lower psi pressure get triggered easily, often without a legitimate reason (e.g., when a tubing extension is added for patient's mobility). Alarm at 6 psi goes off only when a legitimate occlusion occurs.

ORDERING INFORMATION

DESCRIPTION	MFR #	QTY
Bubble Humidifier, 6 psi	32630	1ea/pk, 50pk/cs

Satisfaction Guaranteed

If you are not completely satisfied with any McKesson Brands product, you may return it for a full refund or credit.

Never Compromise on CareSM

With McKesson Brands, you'll never have to compromise. We deliver quality products at competitive prices, backed by our 100% satisfaction guarantee. We have decades of experience in the market, and we offer you support and work to help ensure you have the products you need, when you need them—making it even easier for you to provide the best care for your patients.

The product information contained in this document, including the product images and additional product materials, was collected from various supplier sources. All product claims and specifications are those of the product suppliers and have not been independently verified by McKesson Medical-Surgical or its affiliates ("McKesson"). McKesson is not responsible for errors or omissions in the product information.

The properties of a product may change or be inaccurate following the posting or printing of the product information in the document, either in the print or online version. Caution should be exercised when using or purchasing any products from McKesson's online or print documents by closely examining the product packaging and the labeling prior to use.

Due to product changes, information listed in this document is subject to change without notice. This information is placed solely for your convenience in ordering and McKesson disclaims all responsibility for its completeness and accuracy, whether or not the inaccuracy or incompleteness is due to fault or error by McKesson.