

Super Absorbent Dressings

Super Absorbent Dressings: Super Absorbent Polymers bind drainage and convert it to a gel. These dressings hold more fluid, lock fluid inside the dressing, helping to decrease the risk of maceration and keep the harmful wound fluid components away from the wound and peri-wound skin. Contains protective outer layer and non-adherent wound contact layer for easy removal. Can be used as a primary or secondary dressing for moderate to heavy exudate.

Indications: Primary dressing or secondary dressing for partial thickness or full thickness wounds. Use when there is moderate to large amounts of drainage and extra absorbency is needed to help better manage the fluid and decrease the frequency of dressing changes. Excellent dressing for use under compression for heavily draining venous leg ulcers.

Contraindications: Patients with a known hypersensitivity to the product itself or to its components. Not for use in tunneling wound pockets as the dressing can expand considerably when wound exudates are absorbed. Not for use on wounds with little to no exudates, as this may cause drying out of the wound.

Change Frequency: Change dressing as directed by a healthcare professional.

Recommended Secondary Dressings: Elastic net retainer dressing or conforming bandage. Medical adhesive tape should not be used.



ORDERING INFORMATION

DESCRIPTION	MFR #	QTY	HCPCS	EXUDATE LEVEL*
Super Absorbent Dressing, Sterile, 4" x 5"	61-89545	1/ea, 10/bx, 10bx/cs	A6252	●●●
Super Absorbent Dressing, Sterile, 6" x 9"	61-89569	1/ea, 10/bx, 5bx/cs	A6253	●●●

EXUDATE LEVEL KEY: ●○○ = Light ●●○ = Moderate ●●● = Heavy

*Exudate level indicates highest level of exudate that the dressing can accommodate

Satisfaction Guaranteed
 If you are not completely satisfied with any McKesson Brands product, you may return it for a full refund or credit.

Never Compromise on CareSM

With McKesson Brands, you'll never have to compromise. We deliver quality products at competitive prices, backed by our 100% satisfaction guarantee. We have decades of experience in the market, and we offer you support and work to help ensure you have the products you need, when you need them – making it even easier for you to provide the best care for your patients.

The product information contained in this document, including the product images and additional product materials, was collected from various supplier sources. All product claims and specifications are those of the product suppliers and have not been independently verified by McKesson Medical-Surgical or its affiliates ("McKesson"). McKesson is not responsible for errors or omissions in the product information.

The properties of a product may change or be inaccurate following the posting or printing of the product information in the document, either in the print or online version. Caution should be exercised when using or purchasing any products from McKesson's online or print documents by closely examining the product packaging and the labeling prior to use.

Due to product changes, information listed in this document is subject to change without notice. This information is placed solely for your convenience in ordering and McKesson disclaims all responsibility for its completeness and accuracy, whether or not the inaccuracy or incompleteness is due to fault or error by McKesson.