Gauze Bandage Roll

McKesson Gauze Bandage Rolls are made of 100% prewashed, fluff-dried woven cotton to provide maximum absorbency. Its low-linting formulation helps minimize the risk of fibers being left behind in the wound bed. It may be used as a primary or secondary dressing. Features an open-weave design for fast wicking action, superior aeration and excellent absorbency. Bulk cushioning and conformability improve wound site protection. Ideal for bandaging heads, limbs and difficult-to-dress wounds.

Features and Benefits

- Crinkle-weave pattern increases bulk to cushion and protect wounds
- 6-ply construction with finished edges reduce risk of wound contamination and granulomas
- Open-weave design provides fast wicking action and superior aeration to reduce risk of maceration
- Bulk cushioning and conformability improve wound site protection



ORDERING INFORMATION

DESCRIPTION	MFR #	QТY
Gauze Bandage Roll, Sterile, 2.25" x 3 yds.	16-4062	96 pk/cs
Gauze Bandage Roll, Non-Sterile, 3.4" x 3.6 yds.	16-4063	96 pk/cs
Gauze Bandage Roll, Non-Sterile, 4.5" x 4.1 yds.	30642000	100/cs
Gauze Bandage Roll, Non-Sterile, 2.25" x 3 yds.	16-4262	12/bg, 8 bgs/cs
Gauze Bandage Roll, Sterile, 3.4" x 3.6 yds.	16-4263	1 pk, 96 pk/cs
Gauze Bandage Roll, Sterile, 4.5" x 4.1 yds.	16-4264	1 pk, 100 pk/cs

Satisfaction Guaranteed

If you are not completely satisfied with any McKesson Brands product, you may return it for a full refund or credit.

Never Compromise on Care SM

With McKesson Brands, you'll never have to compromise. We deliver quality products at competitive prices, backed by our 100% satisfaction guarantee. We have decades of experience in the market, and we offer you support and work to help ensure you have the products you need, when you need them — making it even easier for you to provide the best care for your patients.

The product information contained in this document, including the product images and additional product materials, was collected from various supplier sources. All product claims and specifications are those of the product suppliers and have not been independently verified by McKesson Medical-Surgical or its affiliates ("McKesson"). McKesson is not responsible for errors or omissions in the product information.

The properties of a product may change or be inaccurate following the posting or printing of the product information in the document, either in the print or online version. Caution should be exercised when using or purchasing any products from McKesson's online or print documents by closely examining the product packaging and the labeling prior to use.

Due to product changes, information listed in this document is subject to change without notice. This information is placed solely for your convenience in ordering and McKesson disclaims all responsibility for its completeness and accuracy, whether or not the inaccuracy or incompleteness is due to fault or error by McKesson.

© 2017 McKesson Medical-Surgical Inc.