

Hydrocolloid Dressings

Hydrocolloids: Primary or secondary dressing for wounds with light to moderate exudate. Helps maintain a moist wound environment to assist in autolytic debridement. When in contact with wound exudate, it forms a fluid/gel environment over the wound bed. Dressing does not break down or leave residue in the wound.

Thin and Film Back versions feature a smooth, tapered design to help minimize friction and prevent dressing edges from rolling up. Thin hydrocolloid is semi-transparent, allowing visualization of the wound. Contains a semi-transparent layer that permits wound evaluation without removing dressing. Moisture-resistant backing helps prevent urine, feces or moisture from penetrating the dressing to reduce risk of contamination.

Indications: Management of partial and full thickness wounds such as leg ulcers and pressure sores. Help maintain a moist wound healing environment to assist in autolytic debridement and protect the wound from external contamination.

Contraindications: Patients with known hypersensitivities to the components of the dressing or third-degree burns. Conditions that lead to the formation of the wound (pressure, venous insufficiency, etc.) must be considered and addressed as part of the total wound care treatment plan.

Change Frequency: Dressing change frequency is dependent upon the amount of drainage. May be left in place for up to seven days or as directed by a healthcare professional.

Recommended Secondary Dressings: n/a



Satisfaction Guaranteed
 If you are not completely satisfied with any McKesson Brands product, you may return it for a full refund or credit.

Never Compromise on CareSM

With McKesson Brands, you'll never have to compromise. We deliver quality products at competitive prices, backed by our 100% satisfaction guarantee. We have decades of experience in the market, and we offer you support and work to help ensure you have the products you need, when you need them – making it even easier for you to provide the best care for your patients.

The product information contained in this document, including the product images and additional product materials, was collected from various supplier sources. All product claims and specifications are those of the product suppliers and have not been independently verified by McKesson Medical-Surgical or its affiliates ("McKesson"). McKesson is not responsible for errors or omissions in the product information.

The properties of a product may change or be inaccurate following the posting or printing of the product information in the document, either in the print or online version. Caution should be exercised when using or purchasing any products from McKesson's online or print documents by closely examining the product packaging and the labeling prior to use.

Due to product changes, information listed in this document is subject to change without notice. This information is placed solely for your convenience in ordering and McKesson disclaims all responsibility for its completeness and accuracy, whether or not the inaccuracy or incompleteness is due to fault or error by McKesson.

ORDERING INFORMATION

DESCRIPTION	MFR #	QTY	HCPCS	EXUDATE LEVEL*
Hydrocolloid Dressing, Thin, Sterile, 2" x 2"	1882	1/ea, 20/bx	A6234	●○○
Hydrocolloid Dressing, Thin, Sterile, 4" x 4"	1883	1/ea, 10/bx	A6234	●○○
Hydrocolloid Dressing, Thin, Sterile, 6" x 6"	1884	1/ea, 10/bx	A6235	●○○
Hydrocolloid Dressing, Thin, Sacral, Sterile, 6" x 7"	1885	1/ea	A6236	●○○
Hydrocolloid Dressing, Film Backing, Sterile, 2" x 2"	1886	1/ea, 20/bx	A6234	●●○
Hydrocolloid Dressing, Film Backing, Sterile, 4" x 4"	1887	1/ea, 10/bx	A6234	●●○
Hydrocolloid Dressing, Film Backing, Sacral, Sterile, 6" x 7"	1888	1/ea, 5/bx	A6235	●●○
Hydrocolloid Dressing, Foam Backing, Sterile, 4" x 4"	1889	1/ea, 10/bx	A6234	●●○
Hydrocolloid Dressing, Foam Backing, Sterile, 6" x 6"	1890	1/ea, 10/bx	A6235	●●○

EXUDATE LEVEL KEY: ●○○ = Light ●●○ = Moderate ●●● = Heavy

*Exudate level indicates highest level of exudate that the dressing can accommodate