Adhesive and Non-Adhesive Foam Dressings

Adhesive and Non-Adhesive Foam Dressings: Primary or secondary dressing for wounds with moderate to heavy exudate. Hydrocellular foam dressings feature a layer that effectively absorbs wound drainage. Helps maintain a moist healing environment at the wound surface.

The adhesive and non-adhesive with film backing dressings feature a waterproof and breathable outer layer that protects the wound from outside contamination and prevents strike-through of exudate.

Indications: Wounds with moderate to heavy exudate.

Contraindications: Patients with a known hypersensitivity to the product itself or to its components.

Change Frequency: Change dressing as directed by a healthcare professional.

Recommended Secondary Dressings:

Adhesive Foams: n/a

Non-Adhesive Foams: Cover with an appropriate secondary dressing depending on the amount of exudate and frequency of dressing change or secure with medical tape.



Satisfaction Guaranteed

If you are not completely satisfied with any McKesson Brands product, you may return it for a full refund or credit.

Never Compromise on Care SM

With McKesson Brands, you'll never have to compromise. We deliver quality products at competitive prices, backed by our 100% satisfaction guarantee. We have decades of experience in the market, and we offer you support and work to help ensure you have the products you need, when you need them — making it even easier for you to provide the best care for your patients.

The product information contained in this document, including the product images and additional product materials, was collected from various supplier sources. All product claims and specifications are those of the product suppliers and have not been independently verified by McKesson Medical-Surgical or its affiliates ("McKesson"). McKesson is not responsible for errors or omissions in the product information.

The properties of a product may change or be inaccurate following the posting or printing of the product information in the document, either in the print or online version. Caution should be exercised when using or purchasing any products from McKesson's online or print documents by closely examining the product packaging and the labeling prior to use.

Due to product changes, information listed in this document is subject to change without notice. This information is placed solely for your convenience in ordering and McKesson disclaims all responsibility for its completeness and accuracy, whether or not the inaccuracy or incompleteness is due to fault or error by McKesson.

© 2020 McKesson Medical-Surgical Inc.

MSKESSON

ORDERING INFORMATION

MFR #	DESCRIPTION	HCPCS	QTY
16-4670	Hydrocellular Foam Dressings, Acrylic Adhesive, Sterile, 3" x 3"	A6212	1/ea, 10/bx, 10bx/cs
16-4671	Hydrocellular Foam Dressings, Acrylic Adhesive, Sterile, 4" x 4"	A6212	1/ea, 10/bx, 10bx/cs
16-4672	Hydrocellular Foam Dressings, Acrylic Adhesive, Sterile, 6" x 6"	A6213	1/ea, 10/bx, 10bx/cs
16-4673	Hydrocellular Foam Dressings, Acrylic Adhesive, Sterile, 7" x 7"	A6213	1/ea, 10/bx, 10bx/cs
16-4737	Hydrocellular Foam Dressings, Non-Adhesive, Sterile, Film Backing, 4" x 4"	A6209	1/ea, 10/bx, 10bx/cs
16-4738	Hydrocellular Foam Dressings, Non-Adhesive, Sterile, Film Backing, 6" x 6"	A6210	1/ea, 10/bx, 10bx/cs
16-4739	Hydrocellular Foam Dressings, Non-Adhesive, Sterile, 4" x 4"	A6209	1/ea, 10/bx, 10bx/cs

McKesson Medical-Surgical 9954 Mayland Drive Richmond, VA 23233 mms.mckesson.com/mckesson-brands

Post-Acute Care and Home Health 800.347.2456

Home Medical Equipment 888.822.8111

Primary Care 866.625.2679

McKesson Brands Support 800.777.4908