

Laparotomy Sponges

McKesson Laparotomy Sponges enable operating room staff to perform abdominal surgical procedures with precision and ease. Their 4-ply gauze quickly absorbs fluids around the surgical field for better visibility and helps surgeons readily apply pressure to stop heavy bleeding. Dual blue barium sulfate strands are attached to the moisture-retentive pads for easy detection under X-ray equipment. Our laparotomy sponges undergo a vacuum de-linting process that ensures problem-free use by minimizing loose threads and lint.

Features and Benefits

- 100% cotton tightly woven into a prewashed, absorbent design with folded edges for less unraveling
- Available in two sizes for desired selection based on the surgical situation
- X-ray detectable, sterile and non-reusable for patient safety
- Not made with natural rubber latex



ORDERING INFORMATION

DESCRIPTION	MFR #	QTY
Lap Sponges, X-ray Detectable, 12" x 12"	16-2112121	5ea/pk, 40pk/cs
Lap Sponges, X-ray Detectable, 18" x 18"	16-2118181	5ea/pk, 40pk/cs
Lap Sponges, X-ray Detectable, 18" x 18"	16-2118182	5ea/tray, 20trays/cs

Satisfaction Guaranteed

If you are not completely satisfied with any McKesson Brands product, you may return it for a full refund or credit.

Never Compromise on CareSM

With McKesson Brands, you'll never have to compromise. We deliver quality products at competitive prices, backed by our 100% satisfaction guarantee. We have decades of experience in the market, and we offer you support and work to help ensure you have the products you need, when you need them—making it even easier for you to provide the best care for your patients.

The product information contained in this document, including the product images and additional product materials, was collected from various supplier sources. All product claims and specifications are those of the product suppliers and have not been independently verified by McKesson Medical-Surgical or its affiliates ("McKesson"). McKesson is not responsible for errors or omissions in the product information.

The properties of a product may change or be inaccurate following the posting or printing of the product information in the document, either in the print or online version. Caution should be exercised when using or purchasing any products from McKesson's online or print documents by closely examining the product packaging and the labeling prior to use.

Due to product changes, information listed in this document is subject to change without notice. This information is placed solely for your convenience in ordering and McKesson disclaims all responsibility for its completeness and accuracy, whether or not the inaccuracy or incompleteness is due to fault or error by McKesson.