
MSKESSON

Contact Layer Dressings

McKesson contact layer dressings are sterile, versatile dressings that can be used on a variety of different wounds when used with an outer absorbent dressing. The thin transparent dressing is made from polyurethane film and silicone gel with an open perforated structure. These dressings conform to the shape of the wound, while the open perforated structure allows exudate to pass vertically into the outer absorbent dressing and enables easy delivery of topical treatments. They prevent the outer absorbent dressing from sticking to the wound or the skin around the wound. They can also be used as protective layer on superficial wounds that are not discharging and on areas of skin that are very fragile.

Features and Benefits

- Channels allow exudate to transfer vertically into a secondary absorbent dressing and aid delivery of topical treatments
- Designed to allow wound fluid to easily transfer into a secondary dressing
- Protects the wound and the skin
- Prevents an outer dressing from sticking to the wound
- Because the dressing is translucent, it allows the caregiver to asses the healing progress without having to remove the dressing
- Not made with natural rubber latex



ORDERING INFORMATION

| MFR # | DESCRIPTION | QTY |
|-------|--|-------------------------|
| 4816 | Contact Layer Dressings One-Sided, Silicone Adhesive, Sterile, 3 in x 4 in | 1 ea, 10 ea/bx, 4 bx/cs |
| 4817 | Contact Layer Dressings One-Sided, Silicone Adhesive, Sterile, 4 in x 7.2 in | 1 ea, 10 ea/bx, 7 bx/cs |

Satisfaction Guaranteed

If you are not completely satisfied with any McKesson Brands product, you may return it for a full refund or credit.

Never Compromise on CareSM

With McKesson Brands, you'll never have to compromise. We deliver quality products at competitive prices, backed by our 100% satisfaction guarantee. We have decades of experience in the market, and we offer you support and work to help ensure you have the products you need, when you need them—making it even easier for you to provide the best care for your patients.

The product information contained in this document, including the product images and additional product materials, was collected from various supplier sources. All product claims and specifications are those of the product suppliers and have not been independently verified by McKesson Medical-Surgical or its affiliates ("McKesson"). McKesson is not responsible for errors or omissions in the product information.

The properties of a product may change or be inaccurate following the posting or printing of the product information in the document, either in the print or online version. Caution should be exercised when using or purchasing any products from McKesson's online or print documents by closely examining the product packaging and the labeling prior to use.

Due to product changes, information listed in this document is subject to change without notice. This information is placed solely for your convenience in ordering and McKesson disclaims all responsibility for its completeness and accuracy, whether or not the inaccuracy or incompleteness is due to fault or error by McKesson.

MSKESSON

McKesson Medical-Surgical 9954 Mayland Drive Richmond, VA 23233 mms.mckesson.com/mckesson-brands

Post-Acute Care and Home Health 800.347.2456

Home Medical Equipment 888.822.8111

Primary Care 866.625.2679

McKesson Brands Support 800.777.4908