

Collagen Dressings

McKesson collagen dressings are topically-applied advanced wound dressings made from collagen, sodium alginate, carboxymethylcellulose (CMC) and ethylenediaminetetraacetic acid (EDTA). Alginate and CMC allow the dressing to absorb wound drainage, while maintaining an optimal moisture balance to aid in the formation of granulation tissue and epithelialization. The EDTA in the dressing removes zinc to inhibit the activity of matrix metalloproteinases (MMP), thus creating a suitable environment for wound healing.

They are indicated for management of full and partial thickness wounds including pressure ulcers, diabetic ulcers, ulcers caused by mixed vascular etiologies, venous ulcers, donor and graft sites, abrasions, traumatic wounds healing by secondary intention, dehisced surgical wounds, first- and second-degree burns.

Features and Benefits

- Collagen and EDTA work together to target and deactivate excess MMPs, thus triggering the healing process in chronic wounds
- Sheet transforms into gel upon contact with exudates
- Cool, soothing effect of gelling action enhances patient comfort
- May be trimmed to exact wound size
- May be layered for the management of deep wounds
- Sterile, individually-wrapped sheets
- Not made with natural rubber latex



ORDERING INFORMATION

MFR #	DESCRIPTION	QTY
16-1850	Collagen Dressings, Sterile, 2 x 2 in	1 ea, 10 ea/bx, 10 bx/cs
16-1852	Collagen Dressings, Sterile, 4 x 4 in	1 ea, 10 ea/bx, 10 bx/cs

Satisfaction Guaranteed

If you are not completely satisfied with any McKesson Brands product, you may return it for a full refund or credit.

Never Compromise on CareSM

With McKesson Brands, you'll never have to compromise. We deliver quality products at competitive prices, backed by our 100% satisfaction guarantee. We have decades of experience in the market, and we offer you support and work to help ensure you have the products you need, when you need them—making it even easier for you to provide the best care for your patients.

The product information contained in this document, including the product images and additional product materials, was collected from various supplier sources. All product claims and specifications are those of the product suppliers and have not been independently verified by McKesson Medical-Surgical or its affiliates ("McKesson"). McKesson is not responsible for errors or omissions in the product information.

The properties of a product may change or be inaccurate following the posting or printing of the product information in the document, either in the print or online version. Caution should be exercised when using or purchasing any products from McKesson's online or print documents by closely examining the product packaging and the labeling prior to use.

Due to product changes, information listed in this document is subject to change without notice. This information is placed solely for your convenience in ordering and McKesson disclaims all responsibility for its completeness and accuracy, whether or not the inaccuracy or incompleteness is due to fault or error by McKesson.

McKesson Medical-Surgical
9954 Mayland Drive
Richmond, VA 23233
mms.mckesson.com/mckesson-brands

Post-Acute Care and Home Health
800.347.2456

Home Medical Equipment
888.822.8111

Primary Care
866.625.2679

McKesson Brands Support
800.777.4908