

X-Ray Detectable Sponge

McKesson X-Ray Detectable Sponges are suitable for operating room and some surgical applications. A heat-bonded barium sulfate marker is woven into the gauze to help facilitate vivid visibility under X-ray equipment. The gauze sponge features folded edges to resist unraveling and help ensure against the possibility of loose threads. This product's barium sulfate content is greater than 55%, which exceeds requirements.

Features and Benefits

- Made from soft absorbent fabric designed to be virtually lint-free
- 100% prewashed cotton helps ensure high absorbency
- Contains blue X-ray detectable barium sulfate strands for easy identification
- Folded edges to help prevent unraveling and loose threads



ORDERING INFORMATION

DESCRIPTION	MFR #	QTY
X-Ray Detectable Sponge, Sterile, 4" x 4", 16-Ply	16-42446	10/bx, 72bx/cs
X-Ray Detectable Sponge, Sterile, 4" x 8", 16-Ply	16-42447	10/bx, 80bx/cs
X-Ray Detectable Sponge, Sterile, 4" x 8", 12-Ply	16-42846	10/bx, 48bx/cs

Satisfaction Guaranteed

If you are not completely satisfied with any McKesson Brands product, you may return it for a full refund or credit.

Never Compromise on CareSM

With McKesson Brands, you'll never have to compromise. We deliver quality products at competitive prices, backed by our 100% satisfaction guarantee. We have decades of experience in the market, and we offer you support and work to help ensure you have the products you need, when you need them – making it even easier for you to provide the best care for your patients.

The product information contained in this document, including the product images and additional product materials, was collected from various supplier sources. All product claims and specifications are those of the product suppliers and have not been independently verified by McKesson Medical-Surgical or its affiliates ("McKesson"). McKesson is not responsible for errors or omissions in the product information.

The properties of a product may change or be inaccurate following the posting or printing of the product information in the document, either in the print or online version. Caution should be exercised when using or purchasing any products from McKesson's online or print documents by closely examining the product packaging and the labeling prior to use.

Due to product changes, information listed in this document is subject to change without notice. This information is placed solely for your convenience in ordering and McKesson disclaims all responsibility for its completeness and accuracy, whether or not the inaccuracy or incompleteness is due to fault or error by McKesson.