

SmartLink® Desktop 3.0

HCPCS:



- SmartLink® Patient Therapy Management System provides an integrated solution for managing your IntelliPAP® CPAP patients and compliance information, while helping to improve workplace efficiency
- SmartLink Dashboard a collection of customizable tables to quickly identify important tasks, information and patients that need closer attention
- Advanced Find provides power search function to locate patients, doctors,
- payors and devices in the system, and users can save frequently used
- search queries
- Customizable Reports Report templates can be created and saved for reports viewed exactly as you want them, in the order you want them
- Built-in SmartCode Reports SmartCode reports, a DeVilbiss exclusive, are now part of the SmartLink software SmartCode compliance data is downloaded in the SD card with SmartLink and Oximetry data, and can even be viewed over time for people who need more than 90 days to ensure compliance
- Oximetry Reporting Detailed oximetry reports included with the patient report section When using the optional Nonin oximeter, scoring of
- desaturations and heart rate is now part of the patient record

PRODUCTS	Description
SL3-USB	SmartLink 3.0 USB

SPECIFICATIONS	Description
Hard Disk	Up to 2 GB of available space is recommended
Monitor Display	1024 x 768 high color, 32-bit (recommended)
Operating System	Microsoft Windows 7, Windows Vista or Windows 10
Processor	400 MHz Pentium processor or equivalent (minimum)
RAM	512 MB (minimum), 1 GB (recommended)
Software	Adobe Reader, Release 6.0 or greater
Warranty	3 years on module

ACCESSORIES	Description
DV5FC-10	SmartLink Data Card (PK/10), PKG
DV51D-694	SMARTLINK USB EXTENSION CABLE, EA